

Store at 4°C

Annexin V Binding Buffer



天津三箭生物技术股份有限公司
Tianjin Sungene Biotech Co., Ltd.
精准 高效 稳定 Precision Efficient Stable

Catalog Number	Vial Size
AB2000-G	25 tests
AB2000-H	100 tests

Market | 400-621-0003
marketing@sungenebiotech.com

Support | 022-66211636-8024
techsupport@sungenebiotech.com

Web | www.sungenebiotech.com

Important Note: Each Buffer should be diluted by using the same pH PBS as mentioned below.
This product is guaranteed up to one year from purchase.

Product Information

Storage	Keep as concentrated solution. Store at 4°C
Application	FC
Dilution	PBS pH 7.4
Usage	Annexin V Binding Buffer is recommended for use with Annexin V reagents. Cells should be resuspended in Annexin V Binding Buffer at a concentration 1×10^6 /ml and then combined with Annexin V reagents at the recommended volumes.

For Research Use Only.