

Biotin Rat IgG2a Isotype Antibody



天津三箭生物技术股份有限公司
Tianjin Sungene Biotech Co., Ltd.
精准 高效 稳定 Precision Efficient Stable

Catalog Number	Vial Size
R20021-08B	50 µg
R20021-08D	200 µg

Market | 400-621-0003
marketing@sungenebiotech.com

Support | 022-66211636-8024
techsupport@sungenebiotech.com

Web | www.sungenebiotech.com

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.
This product is guaranteed up to one year from purchase.

Purified Antibody Characterization

Clone	Isotype
2A3	Rat IgG2a

Description

Rat IgG2a was chosen as an isotype control after screening on a variety of resting, activated, live, and fixed mouse, rat and human tissues.

Product Information

Conjugation: Biotin

Formulation: PBS pH 7.2, 0.09% NaN₃, 0.2% BSA

Storage: Keep as concentrated solution. Store at 4°C and protected from prolonged exposure to light. **Do not freeze.**

Usage: Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis as negative control. Use at concentrations comparable to those of the specific antibody of interest.

Application: Recommended Application: FC

References

- [1] Barclay, A.N., et al. 1997. The Leukocyte Antigen FactsBook, Academic Press.
- [2] Kishimoto, T.K., et al. 1990. P. Natl. Acad. Sci. USA 87:2244.
- [3] Tedder, T.F., et al. 1995. J. Exp. Med. 181:2259.

For Research Use Only.