

Store at -20°C

## Goat anti-Mouse IgG(H+L)-HRP

Catalog Number

LK2003



天津三箭生物技术股份有限公司  
Tianjin Sungene Biotech Co., Ltd.

精准 高效 稳定 Precision Efficient Stable

**Market** | 400-621-0003  
marketing@sungenebiotech.com

**Support** | 022-66211636-8024  
techsupport@sungenebiotech.com

**Web** | www.sungenebiotech.com

**Important Note:** Centrifuge before opening to ensure complete recovery of vial contents.  
This product is guaranteed up to one year from purchase.

### Purified Antibody Characterization

**Applications:** The source of resistance for an anti-immune mice, Western blot, ELISA and other experiments.

### Description

The secondary antibody is improved based on the original production and have a strong reaction with the IgG heavy and light chains. Meanwhile, the secondary antibody has low reaction with human IgG through adsorption human serum protein. Compared with ordinary secondary antibody, the advantages of this antibody are high sensitivity and specificity and so on. This antibody can apply in western blot (WB), enzyme-linked immunosorbent assay (ELISA) and other experiments.

### Product Information

**Recommended dilutions:** WB: 1/3000-1/5000  
ELISA: 1/10000

**Size(1mg/ml):** 100 $\mu$ l / 1ml

**Storage:** Store at -20°C or below. Avoid repeated freezing and thawing cycles.

**For Research Use Only.**