

Store at -20°C

Goat anti-Rabbit IgG(H+L)-HRP

Catalog Number

LK2001



天津三箭生物技术股份有限公司
Tianjin Sungene Biotech Co., Ltd.

精准 高效 稳定 Precision Efficient Stable

Market | 400-621-0003
marketing@sungenebiotech.com

Support | 022-66211636-8024
techsupport@sungenebiotech.com

Web | www.sungenebiotech.com

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.
This product is guaranteed up to one year from purchase.

Purified Antibody Characterization

Applications: The source of resistance for an anti-immune rabbit, Western blot, ELISA and other experiments.

Description

The secondary antibody is improved based on the original production and have a strong reaction with the IgG heavy and light chains. Meanwhile, the secondary antibody has low reaction with human IgG through adsorption human serum protein. Compared with ordinary secondary antibody, the advantages of this antibody are high sensitivity and specificity and so on. This antibody can apply in western blot (WB), enzyme-linked immunosorbent assay (ELISA) and other experiments.

Product Information

Recommended dilutions: WB: 1/3000-1/5000
ELISA: 1/10000

Size(1mg/ml): 100 μ l / 1ml

Storage: Store at -20°C or below. Avoid repeated freezing and thawing cycles.

For Research Use Only.