

IF488 Goat Anti-Rat IgG(H+L)

| Catalog Number | Vial Size |
|----------------|-----------|
| GR201G-38C | 100 ug |
| GR201G-38E | 500 ug |



天津三箭生物技术股份有限公司
Tianjin Sungene Biotech Co., Ltd.
标准 高效 稳定 Precision Efficient Stable

Market | 400-621-0003
marketing@sungenebiotech.com

Support | 022-66211636-8024
techsupport@sungenebiotech.com

Web | www.sungenebiotech.com

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.
This product is guaranteed up to one year from purchase.

Purified Antibody Characterization

| Isotype | Reactivity |
|---------------------|------------|
| Polyclonal Goat IgG | Rabbit |

Description

Preparation: Purified from goat antiserum by rat IgG affinity chromatography.

Purity: >90% by SDS-PAGE

Specificity: The antibody reacts with whole molecule Rat IgG. No antibody was detected against non-immunoglobulin serum proteins. The antibody may cross-react with immunoglobulins from other species.

Product Information

Conjugation: IF488

Formulation: PBS pH 7.2, 0.09% NaN₃, 0.2% BSA

Concentration: 0.5 mg/ml

Storage: Keep as concentrated solution. Store at 4°C and protected from prolonged exposure to light. **Do not freeze.**

Application: IF(Immunofluorescent)

Usage: IF:1/10-1/1000

The appropriate dilutions must be determined for individual use.

For Research Use Only.