

Store at 4°C

IF647 Goat Anti-Rat IgG(H+L)

Catalog Number	Vial Size
GR201G-40C	100 ug
GR201G-40E	500 ug



天津三箭生物技术股份有限公司
Tianjin Sungene Biotech Co., Ltd.
精准 高效 稳定 Precision Efficient Stable

Market | 400-621-0003
marketing@sungenebiotech.com

Support | 022-66211636-8024
techsupport@sungenebiotech.com

Web | www.sungenebiotech.com

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.
This product is guaranteed up to one year from purchase.

Purified Antibody Characterization

Isotype	Reactivity
Polyclonal Goat IgG	Rabbit

Description

Preparation: Purified from goat antiserum by rat IgG affinity chromatography.

Purity: >90% by SDS-PAGE

Specificity: The antibody reacts with whole molecule Rat IgG. No antibody was detected against non-immunoglobulin serum proteins. The antibody may cross-react with immunoglobulins from other species.

Product Information

Conjugation: IF647

Formulation: PBS pH 7.2, 0.09% NaN₃, 0.2% BSA

Concentration: 0.5 mg/ml

Storage: Keep as concentrated solution. Store at 4°C and protected from prolonged exposure to light. **Do not freeze.**

Application: IF(Immunofluorescent)

Usage: IF:1/10-1/1000

The appropriate dilutions must be determined for individual use.

For Research Use Only.