

Store at 4°C

IF680 Goat Anti-Rabbit IgG(H+L)



天津三箭生物技术股份有限公司
Tianjin Sungene Biotech Co., Ltd.
精准 高效 稳定 Precision Efficient Stable

Catalog Number	Vial Size
GR200G-41C	100 ug
GR200G-41E	500 ug

Market | 400-621-0003
marketing@sungenebiotech.com

Support | 022-66211636-8024
techsupport@sungenebiotech.com

Web | www.sungenebiotech.com

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.
This product is guaranteed up to one year from purchase.

Purified Antibody Characterization

Isotype	Reactivity
Polyclonal Goat IgG	Rabbit

Description

Preparation: The polyclonal antibody was purified from goat antiserum by rabbit IgG affinity chromatography.

Purity: >90% by SDS-PAGE

Specificity: The antibody reacts with whole molecule Rabbit IgG. No antibody was detected against non-immunoglobulin serum proteins. The antibody may cross-react with immunoglobulins from other species.

Product Information

Conjugation: IF680

Formulation: PBS pH 7.2, 0.09% NaN₃, 0.2% BSA

Concentration: 1.0 mg/ml

Storage: Keep as concentrated solution. Store at 4°C and protected from prolonged exposure to light. **Do not freeze.**

Application: IF(Immunofluorescent)

Usage: IF:1/10-1/1000

The appropriate dilutions must be determined for individual use.

iFluor™680 (IF680) has fluorescence excitation and emission maxima close to 682 nm and 701 nm. These spectral characteristics make it excellent alternatives to Cy5.5® and Alexa Fluor® 680 labeling dyes.

For Research Use Only.