

LE/AF Purified Hamster IgG Isotype Antibody

天津三箭生物技术股份有限公司
Tianjin Sungene Biotech Co., Ltd.
精准 高效 稳定 Precision Efficient Stable

Catalog Number	Vial Size
Ha20022-14B	50 µg
Ha20022-14E	500 µg
Ha20022-14F	1 mg

Market | 400-621-0003
marketing@sungenebiotech.com

Support | 022-66211636-8024
techsupport@sungenebiotech.com

Web | www.sungenebiotech.com

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.
This product is guaranteed up to one year from purchase.

Purified Antibody Characterization

Clone	Isotype
Polyclonal	Hamster IgG

Description

The isotype of this antibody is Hamster IgG. This hamster IgG fraction was purified by affinity chromatography from Armenian hamster serum.

Product Information

Production Method: Purification from Armenian Hamster Serum

Medium: n/a

Purification Method: Protein A

Concentration: 1 mg/ml

Endotoxin: < 2.00 EU/mg (LAL)

Purity: > 95%(by SDS-PAGE)

Sterile: 0.2 µm filtration

Formulation: PBS, pH7.2

Storage: Keep as concentrated solution. Store at 4°C as an undiluted liquid. For extended storage aliquot contents and freeze at -20°C or lower. Avoid cycles of freezing and thawing.

For Research Use Only.