

**LE/AF Purified Anti-Mouse CD119 Monoclonal Antibody**

天津三箭生物技术股份有限公司  
Tianjin Sungene Biotech Co., Ltd.  
精准 高效 稳定 Precision Efficient Stable

Catalog Number	Vial Size
M11192-14B	50 µg
M11192-14E	500 µg
M11192-14F	1 mg

<b>Market</b>	400-621-0003 marketing@sungenebiotech.com
<b>Support</b>	022-66211636-8024 techsupport@sungenebiotech.com
<b>Web</b>	www.sungenebiotech.com

**Important Note:** Centrifuge before opening to ensure complete recovery of vial contents.  
This product is guaranteed up to one year from purchase.

**Purified Antibody Characterization**

Clone	Isotype	Reactivity
GR-20	Rat IgG2a	Mouse

**Description**

CDw119 is a 90 kD immunoglobulin superfamily member, also known as IFN- $\gamma$ R $\alpha$  chain. It is a class II cytokine receptor family member that serves as a IFN- $\gamma$ -binding chain associated with the IFN- $\gamma$   $\beta$  chain also known as AF-1. In addition to ligand binding, CDw119 participates in ligand trafficking. CDw119 is expressed on T and B cells, NK cells, fibroblasts, endothelial, and epithelial cells. Binding of IFN- $\gamma$  induces receptor dimerization, internalization, Jak1 and Jak2 protein kinase activation and, ultimately, STAT1 activation. IFN- $\gamma$  initiates and regulates a variety of immune responses.

**Reported Applications**

This GR-20 antibody has been reported for use in flow cytometric staining.

**Product Information**

**Production Method:** Stirred tank fermentation

**Medium:** Hybridoma-SFM + 1%FCS + Gln + Gluc + P/S

**Purification Method:** Protein G

**Concentration:** 1 mg/ml

**Endotoxin:** < 2.00 EU/mg (LAL)

**Purity:** >95% (by SDS-PAGE)

**Sterile:** 0.2 µm Filtration

**Formulated:** PBS, pH7.2

**Storage:** Keep as concentrated solution. Store at 4°C as an undiluted liquid. For extended storage aliquot contents and freeze at -20°C or lower. Avoid cycles of freezing and thawing.

**For Research Use Only.**