

## LE/AF Purified Anti-Mouse CD45R (B220) Monoclonal Antibody



天津三箭生物技术股份有限公司  
Tianjin Sungene Biotech Co., Ltd.  
精准 高效 稳定 Precision Efficient Stable

Catalog Number	Vial Size
MH10456-14B	50 µg
MH10456-14E	500 µg
MH10456-14F	1 mg

**Market** | 400-621-0003  
marketing@sungenebiotech.com

**Support** | 022-66211636-8024  
techsupport@sungenebiotech.com

**Web** | www.sungenebiotech.com

**Important Note:** Centrifuge before opening to ensure complete recovery of vial contents.  
This product is guaranteed up to one year from purchase.

### Purified Antibody Characterization

Clone	Isotype	Reactivity
RA3-6B2	Rat IgG2a	Mouse

### Description

CD45R, also known as B220, is an isoform of CD45. It is a member of the protein tyrosine phosphatase (PTP) family with a molecular weight of approximately 180-240 kD. CD45R is expressed on B cells (at all developmental stages from pro-B cells through mature B cells), activated B cells, and subsets of T and NK cells. CD45R (B220) is also expressed on a subset of abnormal T cells involved in the pathogenesis of systemic autoimmunity in MRL-Fas<sup>lpr</sup> and MRL-Fas<sup>gld</sup> mice. It plays a critical role in TCR and BCR signaling. The primary ligands for CD45 are galectin-1, CD2, CD3, and CD4. CD45R is commonly used as a pan-B cell marker; however, CD19 may be more appropriate for B cell specificity.

### Reported Applications

This RA3-6B2 antibody has been reported for use in Flow Cytometric Analysis, Immunohistochemical and Immunoprecipitation. It has also been reported for use in cell activation.

### Product Information

**Production Method:** Stirred tank fermentation

**Medium:** Hybridoma-SFM + 1%FCS + Gln + Gluc + P/S

**Purification Method:** Protein G

**Concentration:** 1 mg/ml

**Endotoxin:** < 2.00 EU/mg (LAL)

**Purity:** >95% (by SDS-PAGE)

**Sterile:** 0.2 µm Filtration

**Formulated:** PBS, pH7.2

**Storage:** Keep as concentrated solution. Store at 4°C as an undiluted liquid. For extended storage aliquot contents and freeze at -20°C or lower. Avoid cycles of freezing and thawing.

**For Research Use Only.**