

LE/AF Purified Anti-Mouse F4/80 Monoclonal Antibody

天津三箭生物技术股份有限公司
Tianjin Sungene Biotech Co., Ltd.
精准 高效 稳定 Precision Efficient Stable

Catalog Number	Vial Size
M100F1-14B	50 µg
M100F1-14E	500 µg
M100F1-14F	1 mg

Market	400-621-0003 marketing@sungenebiotech.com
Support	022-66211636-8024 techsupport@sungenebiotech.com
Web	www.sungenebiotech.com

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.
This product is guaranteed up to one year from purchase.

Purified Antibody Characterization

Clone	Isotype	Reactivity
BM8	Rat IgG2a	Mouse

Description

The F4/80 monoclonal antibody reacts with mouse F4/80 antigen, an approximately 125 kDa transmembrane protein. The F4/80 antigen is expressed by a majority of mature macrophages and is the best marker for this population of cells. However, other cell types such as Langerhans cells and liver Kupffer cells have been reported to express this antigen. Expression of F4/80 commences during early myeloid development and is upregulated on all BM cells stimulated in vitro with M-CSF. It has been shown that some cytokines downregulate the expression of F4/80 resulting in lack of F4/80 antigen on a subpopulation of macrophages, especially in the lymphoid microenvironment in vivo.

Reported Applications

This BM8 antibody has been reported for use in Flow Cytometric Analysis, Immunohistochemical and Western Blotting.

Product Information

Production Method: Stirred tank fermentation

Medium: Hybridoma-SFM + 1%FCS + Gln + Gluc + P/S

Purification Method: Protein G

Concentration: 1 mg/ml

Endotoxin: < 2.00 EU/mg (LAL)

Purity: >95% (by SDS-PAGE)

Sterile: 0.2 µm Filtration

Formulated: PBS, pH7.2

Storage: Keep as concentrated solution. Store at 4°C as an undiluted liquid. For extended storage aliquot contents and freeze at -20°C or lower. Avoid cycles of freezing and thawing.

For Research Use Only.