

LE/AF Purified Anti-Mouse IL-12/IL-23 p40 Monoclonal Antibody



天津三箭生物技术股份有限公司
Tianjin Sungene Biotech Co., Ltd.
精准 高效 稳定 Precision Efficient Stable

Catalog Number	Vial Size
M100I121-14B	50 µg
M100I121-14E	500 µg
M100I121-14F	1 mg

Market | 400-621-0003
marketing@sungenebiotech.com

Support | 022-66211636-8024
techsupport@sungenebiotech.com

Web | www.sungenebiotech.com

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.
This product is guaranteed up to one year from purchase.

Purified Antibody Characterization

Clone	Isotype	Reactivity
C17.8	Rat IgG2a	Mouse

Description

The C17.8 antibody reacts with mouse IL-12 p40 subunit of the IL-12 p70 heterodimer and IL-23 p40 subunit of the IL-23 p19/p40, as well as p40 monomer and homodimer. The C17.8 antibody can neutralize the bioactivity of natural or recombinant IL-12.

Reported Applications

This C17.8 antibody has been reported for use in Western Blotting, Flow Cytometric Analysis, ELISA (Enzyme-linked Immunosorbent Assay), Immunoprecipitation. It has also been reported for use in Cytokine Neutralization.

Product Information

Production Method: Stirred tank fermentation

Medium: Hybridoma-SFM + 1%FCS + Gln + Gluc + P/S

Purification Method: Protein G

Concentration: 1 mg/ml

Endotoxin: < 2.00 EU/mg (LAL)

Purity: >95% (by SDS-PAGE)

Sterile: 0.2 µm Filtration

Formulated: PBS, pH7.2

Storage: Keep as concentrated solution. Store at 4°C as an undiluted liquid. For extended storage aliquot contents and freeze at -20°C or lower. Avoid cycles of freezing and thawing.

For Research Use Only.