

## PE Goat Anti-Human IgG(H+L)



天津三箭生物技术股份有限公司  
Tianjin Sungene Biotech Co., Ltd.  
精准 高效 稳定 Precision Efficient Stable

Catalog Number	Vial Size
GH200G-09C	100 ug
GH200G-09E	500 ug

**Market** | 400-621-0003  
marketing@sungenebiotech.com

**Support** | 022-66211636-8024  
techsupport@sungenebiotech.com

**Web** | www.sungenebiotech.com

**Important Note:** Centrifuge before opening to ensure complete recovery of vial contents.  
This product is guaranteed up to one year from purchase.

### Purified Antibody Characterization

**Reactivity**  
Human

#### Description

**Preparation:** The Goat Anti-Human IgG(H+L) antibody was purified from goat antiserum by human IgG affinity chromatography.

**Purity:** >90% by SDS-PAGE

**Specificity:** The secondary antibody is improved based on the original production and have a strong reaction with the IgG heavy and light chains.

### Product Information

**Conjugation:** PE

**Formulation:** PBS pH 7.2, 0.09% NaN<sub>3</sub>, 0.2% BSA

**Concentration:** 0.2 mg/ml

**Storage:** Keep as concentrated solution. Store at 4°C and protected from prolonged exposure to light. **Do not freeze.**

**Usage:** Recommended dilution:

IF: 1/50 ~ 1/800

Since applications vary, the appropriate dilutions must be determined for individual use.

**For Research Use Only.**